

# Hyperhidrosis Treatment

Solves persistent heavy perspiration of hands, feet, underarms and crotch

## METHOD: Iontophoresis

(Via continuous and pulsed current therapy)

Proven effective  
Standard approach in Germany for over 30 years already.  
No surgery, No injections, No side effects.  
Also help treat gram negative foot infection

**KaWe SwiSto3**

05.11110.002

**Tigil-Pasma Device**



## FEATURES

- Simple operation and menu navigation
- Choice of continuous or pulsed current
- Automatic current control, no feeling of electric shock when hands and feet are removed from the water trays
- Treatment of the hands without foot switch possible
- Automatic shut-off, sleep mode
- Lio-io rechargeable battery with charging function, charging plug
- Per change: up to 10, 20 minute treatments; Easy to disinfect

*Helps with heavy perspiration/pasma*



**Underarm**

**Feet**

## How to use:

- Cover electrodes with mesh pad
- Fill tray with approx. 3-4cm of lukewarm water enough to cover hands or feet surface.
- If for underarm use ensure sponge pockets fully soaked in water
- 15-20 minutes at a time intially at 3 times per week (approximately 10-15 treatments) to get desired healing results
- Maintenance treatment on atleast once a week
- Adjust current as needed (but max of 25 mA)
- Recommended switch polarity if thickly calloused

*Consult doctor if you have any of the following: Cardiac arrhythmia, Pregnancy, Pacemakers or metal implants, Metal-containing intrauterine devices (contraceptive coil), Major skin defects that cannot be covered with Vaseline, Insensitivity to painful stimuli*