



Six Drawer Mini Line Emergency Cart with Extended Stabilizing Bumper 5" Wheels, Specialty Accessory Package, with Breakaway Lock #3256B-EMG

Harloff's Mini Line of carts offer a narrow design with reduced footprint. Suitable for hospitals, clinics and doctor offices with both space and budget constraints



Specifications

Construction:

- Painted steel cabinet and drawers

Features:

- Space-saver design fits into smaller areas
- Allows supply storage in areas with space limitations
- Cabinet with 30 vertical inches (76.2 cm) drawer space
- Breakaway lock system with 100 plastic breakaway seals
- Bumper system extends wheel base and stabilizes cart
- 5" ball-bearing full swivel casters – one brake, one directional
- Pre-threaded accessory mounting holes
- Ball-bearing full extension drawer slides
- Removable plastic top
- Latex free
- Twelve year warranty
- Durable hard-baked powder coat finish

Accessory Package Includes:

- 2-prong IV pole (IV pole extends to 76.5") (680401)
- Oxygen tank holder (AL680403)
- Cardiac board and brackets (680406)
- Aluminum push handle / side rail (68161)

Drawer Configuration:

- (3) - 3-inch (7.6 cm)
- (2) - 6-inch (15.2 cm)
- (1) - 9-inch (22.9 cm)

Internal Drawer Dimensions:

- 12.6" W x 16.75" D (32cm x 42.5cm)

Assembled Cart Dimensions:

- 69" H x 24.5" W x 25" D (175.3 cm H x 62.2 cm W x 63.5 cm D) height includes IV Pole – extends to 87.5" / 222.25 cm H
- Cart net weight: 106 pounds (48.5 kg)
- Approximate shipping weight: 126 pounds (57.2 kg)

Sixteen available color choices:

*Antimicrobial Finish Available

Sand, Gloss Beige, Burgundy, Light Blue, Navy, Light Gray, Mauve, Red, Teal, Yellow, Hammer Tone Blue, Hammer Tone Gray, Hammer Tone Green, Hammer Tone Red, Hammer Tone Black and White

Exclusively Distributed by:

WellnessPRO
tomorrow's medical devices