



3 Drawer and Storage Compartment GP Line Emergency Cart Specialty Package #6161

Specifications



Construction:

- All Steel

Features include:

- Stainless steel top rail/push handle
- Removable plastic top
- Tall GP cabinet with stainless steel breakaway locking system
- Breakaway lock seals, pkg-50
- Vinyl-wrapped full wrap-around bumper system
- Five-inch (12.7 cm) urethane single-wheel caster set, full swivel – two with brake
- Storage compartment: 18" H (45.7 cm) with locking door and one adjustable shelf
- Ball-bearing drawer slides
- Single-walled drawer faces
- Cart is delivered fully assembled
- Certified "Latex Free"
- Hard-baked powder coat finish
- Compatible with select accessories and drawer dividers
- Five Year Warranty

Drawer configuration:

- 2 x 3-inch (7.6 cm)
- 1 x 6-inch (15.2 cm)

Drawer dimensions:

- 23" W x 17" D (58.4 cm W x 43.2 cm D)

Twelve available color choices:

- Beige (-BG), Burgundy (-BU), Light Blue (-LB), Navy (-NY), Light Gray (-GY), Mauve (-MV), Red (-RD), Teal (-TL), Yellow (-YL), Hammer Tone Blue (-HTDB), Hammer Tone Gray (-HTGY) and Hammer Tone Green (-HTGN)

Assembled cart dimensions:

- 43" H x 28.5" W x 22.25" D (107 cm H x 71 cm W x 55.9 cm D)
- Cart Net Weight: 125 pounds (56.9 kg)
- Approx. shipping weight: 145 pounds (66 kg)