

HARLOFF

Classic Line, Nine Drawer, Pediatric Resuscitation Cart, Specialty Package #6401PEC



**Color Shown is Beige.
Additional Colors Available**

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- Tall cabinet with separate breakaway locks on each drawer
- Pack of 50 break-away seals
- 30" of vertical drawer space
- Cart is designed to be compatible with popular pediatric color coding systems
- Includes a pack of multi-colored labels for the end user to apply to identify drawers
- Stabilizing pontoon bumpers
- 5" ball-bearing casters, all full swivel - one tracking, one with brake
- Removable plastic top
- Painted Steel top rail
- Pre-drilled threaded holes for side mounted accessories
- Ball-bearing full extension drawer slides
- Latex free
- Durable, hard-baked powder coat finish
- Cart is delivered assembled
- Twelve year warranty on manufactured parts

Accessory Package Includes:

- Cardiac board and brackets (680406)
- Drop shelf (680302)
- Utility hooks (4 ea) (680402)
- IV pole (680401)
- Oxygen tank holder - D/E size cylinders (680403)

Assembled Cart Dimensions:

- 45.75" H x 35" W x 22" D
(116.21 cm H x 89 cm W x 55.9 cm D)
- With IV pole: 86.25" (219 cm H)
- Cart weight: 172 lbs (78 kg)
- Approximate shipping weight: 213 lbs (96.6 kg)
- Approximate shipping dims: 53" H x 33" W x 26" D
(134.6 cm H x 84 cm W x 66 cm D)

Drawer Configuration:

- (8) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)