

HARLOFF

Classic Line, Super-Stat Emergency Cart # 6411



**Color Shown is Red.
Additional Colors Available**

Internal Drawer Dimensions:

23" W x 17" D (58.4 cm W x 43.2 cm D)

Assembled Cart Dimensions:

- 45.75" H x 39" W x 24" D
(116.2 cm H x 99 cm W x 61 cm D)
- With IV pole: 86.25" (219 cm)
- Cart net weight: 181 lbs (82.1 kg)
- Approximate shipping weight: 235 lbs (106.6 kg)
- Approximate shipping dims 53" H x 48" W x 30" D
(134.6 cm H x 121.9 cm W x 76.2 cm D)

Specifications

Construction:

- Sturdy, uni-body steel construction

Features:

- 30" of vertical drawer space
- Painted steel break-away lock bar (30")
- Pack of 100 break-away seals
- Painted steel top rail / push handle
- Ball-bearing full extension drawer slides
- 5" ball-bearing casters; full swivel – one with brake, one tracking
- Removable plastic top
- Stabilizing pontoon bumpers
- Latex free
- Durable, hard-baked powder coat finish
- Twelve year warranty on manufactured parts

Super Stat Accessory Package Includes:

- Cardiac board and brackets (680406)
- Three utility hooks (680402)
- IV pole (680401)
- Painted steel defibrillator shelf 20" x 12" (680323)
- Suction unit shelf (680305)
- Oxygen tank holder - D/E size cylinders (AL680403)
- Chart holder (680417)
- Locking sharps container, 5qt sharps, and glove box holder (684811)
- Two large utility baskets (680447)
- One small utility basket (680431)
- Two 6 compartment drawer trays (fit three per drawer), (40108)
- One 7 compartment drawer tray (fit three per drawer), (40109)
- Two 3" drawer trays with adjustable dividers (fit two per drawer), (681501)
- Four 6" drawer trays with adjustable dividers (fit two per drawer), (681502)
- Hospital Grade electrical Outlet Strip with Cord Wrap not included as of January 15, 2015 due to UL specification change. Outlet conforming to UL1363A can be purchased as an accessory. (680404A)

Drawer Configuration:

- (1) – 6"
- (3) – 3"
- (1) – 6"
- (1) – 9"

Exclusively Distributed by:

WellnessPRO
tomorrow's medical devices