

Conveniently giving all the best for your baby

www.mamivac.com

mamivac®

Breast Pump Brand

MADE IN GERMANY 



Sensitive CFH

REF 286000

- Practical stand for mobile use in clinics
- With rechargeable batteries (optional)
- Integrated dual bottle holder
- Continuously adjustable
- Stimulation phase and expressing phase
- Easy to operate

Sensitive C

REF 284000

- The mamivac® 2-phase technology imitates the baby's sucking
- No suction loss with double pumping
- Gentle, hygienic and silent operation guaranteed
- Compact, maintenance-free and durable device
- Germ free – patented protective diaphragm in the suction kit offers the best possible protection and safety for mother and child

Lactive

REF 287000

- Easy to use
- Expressing comfort due to 2-phase technology
- Visual stimulation with a photo of your baby. This stimulus supports the flow of breast milk.
- Time-saving due to optional double pumping
- Best hygienic standard for mother & child

Easy Manual

REF 281360

- The mamivac® 2-phase technology imitates the baby's sucking
- Adjustment valve for continuously adjustable suction strength
- Soft breastshield with massaging effect
- Hygienic as the breast milk flows straight into the bottle
- With the enclosed freezer bag adapter, milk can be expressed straight into a mamivac® FREEZER BAG

Exclusively Distributed by:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Avenue corner Lourdes St., Ortigas Center Pasig 1605, Philippines

Tel: (02) 636-3580 • Fax: (02) 636-7542
Mobile: (63) 920-947-7185, (63) 922-823-2282

Email: wellnessproinc@gmail.com
jcjuan@wellnessproinc.com
www.wellnessproinc.com

mamivac®

Breast Pump Brand

MADE IN GERMANY **KaWeCo**



Other Products:



Technical Data

	Sensitive CFH	Sensitive C	Lactive
Vacuum (suction strength)	Continuously adjustable from -5 to -33 k Pa (± 15%)	Continuously adjustable from -5 to -33 k Pa (± 15%)	Continuously adjustable from -5 to -33 k Pa (± 15%)
Cycle (speed)	Stimulation phase with 100 cycles/min. Continuously adjustable expressing phase with 20-48 cycles/min.	Stimulation phase with 100 cycles/min. Continuously adjustable expressing phase with 20-48 cycles/min.	Stimulation phase 60 cycles / wmin. expressing phase 30 cycles / min.
Power supply	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC 12 V/DC, max. 1.25 A/DC	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC 12 V/DC, max. 1.25 A/DC	100-240 V/AC, 50/60 Hz, max. 0.62 A/AC 12 V/DC, max. 1.25 A/DC
Dimensions H x D x L[]	24 cm x 21 cm x 22 cm	19 cm x 11.2 cm x 14.5 cm	11.5 cm x 12.5 cm x 19.5 cm
Weight	1.3 kg	0.8 kg	0.6 kg
Warranty	3 years	3 years	1 year
Built-in Technologies	2-Phase Technology Double Pumping Protective Diaphragm	2-Phase Technology Double Pumping Protective Diaphragm	2-Phase Technology Double Pumping Protective Diaphragm



2-PHASE TECHNOLOGY

To make the expressing of breast milk as natural as possible, we have developed the mamivac 2-phase technology. The continuously adjustable phases consist of a stimulation phase and an expressing phase.



DOUBLE PUMPING

The simultaneous expressing of milk from both breasts saves time which will reduce express time by about 50%.



PROTECTIVE DIAPHRAGM

This diaphragm hermetically seals off the breastpump. This avoids cross contamination between users and the device itself.

KaWeCo GmbH

Gerlinger Str. 36-38
71254 Ditzingen
Germany

P h o n e: +49 (0) 7156/17602 - 100

F a x: +49 (0) 7156/17602 - 500

info@kaweco.de
www.kaweco.de

www.KaWeCo.de

Exclusively Distributed by:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Avenue corner Lourdes St., Ortigas Center Pasig 1605, Philippines

Tel: (02) 636-3580 • Fax: (02) 636-7542
Mobile: (63) 920-947-7185, (63) 922-823-2282

Email: wellnessproinc@gmail.com
jcjuan@wellnessproinc.com
www.wellnessproinc.com