

Exclusive Philippine Distributor:

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jjuan@wellnessproinc.com
www.wellnessproinc.com

WellnessPRO
tomorrow's medical devices

FAMILIE | PRODUCT LINE KaWe EUROLIGHT®



KaWe EUROLIGHT® VET C30

- Clic-Verschluss
 - Vakuum-Lampe
 - 2-fache Lupenvergrößerung mit spezieller Brennweite
 - Metall-Batterie-/Ladegriff C
 - Rheostat dimmbar
 - 3 VET-Dauer-Ohrtrichter (▶23)
Ø 4,0 / 5,0 / 7,0 mm
 - in Reißverschlussstasche
 - Batterien: 2 Baby (Typ C) ☒
- click-on connection
 - vacuum bulb
 - two-fold magnification with special focal distance
 - metal battery/charging handle C
 - dimmable rheostat
 - 3 reusable funnels (VET) (▶23)
Ø 4.0 / 5.0 / 7.0 mm
 - in bag with zipper
 - batteries: 2 C batteries ☒



4 **VET C30** | 2,5V | VL REF 01.11135.001 (24855)

Leuchtmittel | illuminant REF 12.75111.013 (28943)



KaWe EUROLIGHT® VET C30 OP

- Clic-Verschluss
 - Vakuum-Lampe
 - 2-fache Lupenvergrößerung mit spezieller Brennweite
 - Metall-Batterie-/Ladegriff C
 - Rheostat dimmbar
 - 3 VET-Dauer-Ohrtrichter (▶23)
Ø 4,0 / 5,0 / 7,0 mm
 - in Reißverschlussstasche
 - speziell für den OP-Bereich
 - Otophoskop-Kopf oben zugänglich für Instrumentarium
 - Batterien: 2 Baby (Typ C) ☒
- click-on connection
 - vacuum bulb
 - two-fold magnification with special focal distance
 - metal battery/charging handle C
 - dimmable rheostat
 - 3 reusable funnels (VET) (▶23)
Ø 4.0 / 5.0 / 7.0 mm
 - in bag with zipper
 - specially designed for operational procedures
 - head accessible for instruments
 - batteries: 2 C batteries ☒



5 **VET C30 OP** | 2,5V | VL REF 01.11136.001

Leuchtmittel | illuminant REF 12.75111.013 (28943)



KaWe EUROLIGHT® F.O.30

- Clic-Verschluss
 - Xenon-Lampe
 - Metall-Batterie-/Ladegriff C
 - Rheostat dimmbar
 - Einweg-Ohrtrichter (▶22)
10 x Ø 2,5 mm
10 x Ø 4,0 mm
 - in Reißverschlussstasche
 - Batterien: 2 Baby (Typ C) ☒
- click-on connection
 - xenon bulb
 - metal battery/charging handle C
 - dimmable rheostat
 - disposable ear funnels (▶22)
10 x Ø 2.5 mm
10 x Ø 4.0 mm
 - in bag with zipper
 - batteries: 2 C batteries ☒



6 **F.O.30** | 2,5V | XL REF 01.11330.001 (24750)

Leuchtmittel | illuminant REF 12.75144.013



KaWe EUROLIGHT® F.O.30 LED

- Clic-Verschluss
 - LED-Lampe
 - Beleuchtungsstärke ca. 8.000 Lux
 - Lebensdauer Lampe ca. 100.000 h
 - Metall-Batteriegriff C
 - Einweg-Ohrtrichter (▶22)
10 x Ø 2,5 mm, 10 x Ø 4,0 mm
 - in Reißverschlussstasche
 - Batterien: 2 Baby (Typ C) ☒
- click-on connection
 - LED bulb
 - illuminance approx. 8,000 Lux
 - life of bulb approx. 100,000 h
 - metal battery handle C
 - disposable ear funnels (▶22)
10 x Ø 2.5 mm, 10 x Ø 4.0 mm
 - in bag with zipper
 - batteries: 2 C batteries ☒



7 **F.O.30 LED** | 2,5V | LED REF 01.11530.002

Leuchtmittel | illuminant REF 12.75153.003

aufladbar in Ladestation KaWe MedCharge® 4000 in Verbindung mit Ladebatterie / rechargeable in KaWe MedCharge® 4000 charging station with rechargeable battery
 Batterien sind nicht im Lieferumfang enthalten! / Batteries are not included!

VL Vakuum-Lampe / vacuum bulb
XL Xenon-Lampe / xenon bulb
LED LED-Lampe / LED bulb