

# PRAXIS-EINRICHTUNG | MEDICAL EQUIPMENT

## Praxisliege Classic | Examination couch Classic

- Kunstleder
  - Polsterstärke: 55 mm
  - abwasch- und desinfizierbar
  - Gesamtgewicht: 32 kg
  - Nenntraglast: 400 kg
  - Bewegte Last: 135 kg
  - verstellbares Rückenteil: Rastschere in 9 Stufen verstellbar bis ca. +34°
  - Fingerschutz durch Polsterdistanz am beweglichen Rückenteil
  - Gestell: pulverbeschichtet, grauweiß
  - inkl. Papierrollenhalter
  - Füße mit verfärbungsfreien Nivellier-Stopfen
  - Lieferung: 2 Packstücke
- synthetic leather
  - thickness of cushion: 55 mm
  - washable, disinfectable
  - overall weight: 32 kg
  - rated load: 400kg
  - moving load: 135kg
  - adjustable backrest: 9 latch settings up to approx.+34°
  - cushion distance on movable backrest protects fingers
  - frame: powder-coated, grey-white
  - incl. paper roll holder
  - legs with non-marking leveling feet
  - Delivery: 2 packaged items



171

## Fußauftritt | Step stool

- aus Edelstahl
- 1- oder 2-stufig erhältlich
- Tragfähigkeit: ca. 200 kg
- made of stainless steel
- available with 1 or 2 steps
- load capacity: approx. 200 kg



1-stufig | 1 step 40 x 25 x 22 cm  
2-stufig | 2 steps 40 x 25 x 34 cm

1-stufig   1 step	REF 09.61100.002
2-stufig   2 steps	REF 09.62100.002

## Drehhocker Classic | Swivel stool Classic

- Kunstleder
  - Sitzhöhe: 450 – 570 mm
  - Polsterstärke: 50 mm
  - Sitzpolster Ø 37 cm
  - Gesamtgewicht: ca. 5,5 kg
  - Tragfähigkeit: max. 110 kg
  - Sicherheitsgasfeder mit Hebelauslösung
  - Laufrollen für Teppichböden geeignet
- synthetic leather
  - seat height: 450 – 570 mm
  - thickness of cushion: 50 mm
  - seat cushion Ø 37 cm
  - overall weight: approx. 5.5 kg
  - load capacity: 110 kg max.
  - safety pneumatic spring with release lever
  - the castors are suitable for carpet



schwarz   black	REF 09.51302.022
blau   blue	REF 09.51302.032
grau   grey	REF 09.51302.062

