



Preconfigured Triple Column Medical Storage Cart with Tambour Door #MS8160-VASC

The MS8160-VASC is a triple column cart with tambour door, designed to store medical supplies related to vascular surgery and diseases of the vascular system, featuring the most innovative modular tray system on the market.

The MS8160-VASC is preconfigured as follows: left column, top to bottom, 2ea 1/2" trays, 2ea 2" trays, 3ea 4" trays, 3ea 8" trays; middle column, top down, 2ea 1/2" trays, 2ea 2" trays, 3ea 4"trays, 3ea 8" trays; right column, top to bottom, 2 ea 1/2" trays, 2ea 2" trays, 3ea 4" trays, 3ea 8" trays.



Assembled cart dimensions:

- 80.5" H x 59.6" W x 27.9" D
(204.5 cm H x 151.4 cm W x 70.8 cm D)
- Cart net weight: 400 pounds (181 kg)
- Approximate shipping weight: 590 pounds (267.2 kg)

Construction:

- Painted steel cabinet, 20-gauge powder-coated steel

Advantages and Benefits:

- Offers 30% more tray or basket storage capacity when compared with leading competitor's product, using 2", 4" and 8" modular trays; or 4" and 7" wired baskets.
- Increased storage capability means fewer carts need to be purchased, thereby reducing capital purchase cost.
- Further optimize space within trays with adjustable dividers that have a unique locking design, called the "clip-on effect," to withstand high loads.
- Unique design of plastic side panels, combined with pull-out-stoppers, allows modular trays to be extended and viewed at an angle without falling out of the cart.
- Catheter shelf can be added featuring three robust dual slide assemblies, each slide contains five hooks.

Standard Features:

- One half inch trays: Total of six trays that are 23.50 x 15.75 x 1/2 inches
- Two inch modular trays: Total of six trays that are 23.50 x 15.75 x 2 inches
- Four inch modular trays: Total of nine trays that are 23.50 x 15.75 x 4 inches
- Eight inch modular trays: Total of nine trays that is 23.50 x 15.75 x 8 inches
- Plastic side panels form a modular system consisting of 141 insertion opportunities (47 per column) to hold modular exchange trays and baskets, which are sold separately.
- Catheter Shelf Holders can be added to any column.
- One of the columns can be removed to increase catheter storage capacity
- Roll-up tambour door with key lock allows access to supplies during clinical hours.
- Integrated push handles are flush with cart side, positioned vertically for better control
- Corner bumpers with premium five inch casters, full swivel – two with brake
- Standard color is Sand. Also available in 16 other colors, using durable hard-baked powder coat finish
- Ships fully assembled with twelve year warranty