

# seca 725

## Mechanical baby scales with sliding weights



- Classic design and high-quality precision mechanics.
- Precise ascertainment of weight due to zero adjustment.
- Softly curved weighing tray for maximum safety.
- Easy-care and mobile.

## seca 725:

### The sturdy classic.

**A** baby gains weight quickly – its weight has almost doubled after four to six months. As soon as the baby begins to crawl and run, the one or other gram of baby fat will quickly disappear. In order that the baby's development can be monitored safely from the very beginning, instruments are required in pediatrics on which one can rely long-term. A baby scale that will work in all surroundings: The mechanical baby scale seca 725 has been around for over 50 years and is still there being used by health professionals all over the world.



#### Classic design and high-quality precision mechanics.

The seca 725 combines the classic flat construction with high-quality precision mechanics. With its smooth running weights it is easy to use and provides precise results. The baby lies safe and secure in the softly curved weighing tray. The robust steel construction and cast levers ensure years of durability and accuracy. All this has made this mechanical baby scale a favorite for many decades.

#### Precise weight ascertainment due to zero adjustment.

The zero adjustment ensures the reliable determination of the net weight – a soft mat can simply be tared away, showing the actual weight of the baby. The baby scale has lbs and kg measurements on the balance beam and the weight graduations are very small, 5 g. The weight capacity is 16 kg respectively. These graduations are hard to beat on any mechanical baby scale you can find and because this unit is mechanical, there is very little that can go wrong with it.

#### Easy-care and mobile.

The weighing tray is easy to remove from the base for cleaning. The hard-wearing, scratch-proof finish and the smooth, seamless surface makes the seca 725 extremely easy to clean. The scale can be safely stored and is easy to transport in the optional carrying case seca 425.



To simplify cleaning, weighing tray can be removed.



The reset-to-zero function guarantees reliable determination of net weight by simply deducting the extra weight of padding or diaper.

With plastic wheels and retractable handle.



## seca 725

### Technical Data

- Capacity: 16 kg, 32 lbs
- Graduation: 5 g, 1/4 oz
- Dimensions (WxHxD):  
550 x 180 x 290 mm /  
21.7 x 7.1 x 11.4"
- Dimensions/tray (WxHxD):  
550 x 110 x 270 mm /  
21.7 x 4.3 x 10.6"
- Weight: 6.4 kg / 14.1 lbs
- Functions: reset-to-zero function
- Optional: carrying case seca 425

Exclusive Philippine Distributor:

**WellnessPRO**  
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,  
Ortigas Ctr. Pasig City 1605 Philippines  
Mobile no.: +63 922 896 0964  
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473  
Fax: +63 2 636 7542  
Email: wellnessproinc@gmail.com, icjuan@wellnessproinc.com  
[www.wellnessproinc.com](http://www.wellnessproinc.com)