

seca 952

Chair scale for weighing while seated



- The simple solution for weighing wheelchair patients.
- Fold-up armrests and footrests make sitting down easy.
- Brakes on wheels for secure stance.
- Intuitive operating elements simplify use.
- All essential functions for maximum usage.

seca 952:

Relief for patient and medical personnel.

Standing on a scale is a problem for physically disabled or long-term care patients. Chair scales are of great assistance to patients and medical staff in nursing homes, hospitals, geriatric and orthopedic departments. The seca 952 offers a sound, user-oriented solution with all the functions that a chair scale should offer.

Everything that counts.

With a capacity of 440 lbs and a graduation of 0.2 lbs, the seca 952 fulfills all basic medical standards. Added to that are the typical seca advantages such as fold-up armrests and footrests, secure wheel brakes, an ergonomically shaped seat and smooth running rollers, which, together with the ergonomically shaped handle, make it very easy to push the chair scale.

Simple operation.

Thanks to the fold-up armrests and footrests, elderly and heavy patients can take a seat on the scale easily and comfortably. It is just as easy for personnel to use the scale because the batteries eliminate the need for connection to an electrical outlet and keep on working through up to 4000 weighings. No time is lost to a learning curve – a user can operate the scale intuitively. Display elements are conveniently located in the handle and can be set to show kilograms, pounds or stones.

TARE and auto-HOLD.

Of course the seca 952 is also equipped with the important TARE function which allows every newly added weight to be ascertained separately and the auto-HOLD function which lets medical personnel take care of the patient first and read out the results later. All these advantages plus well-known workmanship and durability make the scale a genuine seca product. It fulfills all the expectations people have of the world leader in medical weighing and measuring.

The new ergonomic handle makes scale easier to push.



To make sitting down and standing up easier, the armrests and footrests fold up out of the way.



The practically positioned display is easy to use.

seca 952

Technical data

- Capacity: 440 lbs / 200 kg
- Graduation: 0.2 lbs / 100 g
- Dimensions (WxHxD):
22.17 x 35.31 x 37.64" / 563 x 897 x 956 mm
- Weight: 49.42 lbs / 22.4 kg
- Power supply: Power adapter
- Functions: TARE, auto-HOLD, lbs/kg switch-over, automatic switch-off

Exclusive Philippine Distributor:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jcjuan@wellnessproinc.com
www.wellnessproinc.com

Toll free +1 800 542 7322 • Fax +1 888 705 7397 • info.us@seca.com • www.seca.com
North America West • seca corp. • 13601 Benson Avenue • Chino, CA 91710 USA
North America East • seca corp. • 7240 Parkway Drive, Suite 120 • Hanover, MD 21076 USA