

HARLOFF

V-Series, Six Drawer Crash Cart, Specialty Accessory Package #V30-6EMG

This high capacity crash cart is designed for healthcare facilities under budget pressure requiring a durable, fully equipped emergency cart offering Harloff's quality construction at a value price.

This emergency cart includes all commonly required accessories and is perfect for hospitals, clinics, dialysis centers, doctors' offices, urgent care facilities and nursing homes needing to buy carts to store, transport and dispense critical supplies used for emergency situations.



Drawer Configuration:

- (4) – 3" (7.6 cm)
- (1) – 6" (15.2 cm)
- (1) – 12" (30.5 cm)

Integrated Pull Handle on all drawers

Sixteen available color choices:

*Antimicrobial Finish Available

Burgundy, Hammer Red, Red, Yellow, Beige, Sand, White, HammerTone Green, Teal, HammerTone Blue, Navy, Light Blue, HammerTone Gray, Light Gray, Mauve, Hammer Black

Specifications

Construction:

- Steel construction

Features:

- Break-away lock bar
- 100 break-away seals included
- 30" of vertical drawer space
- Stabilizing bumper frame with extended front casters and (4) corner bumpers
- 5" (12.7 cm) ball-bearing casters, all full swivel - two with brake
- Pre-drilled threaded holes for V-Series side mounted accessories
- Ball-bearing full extension drawer slides
- Twelve year manufacturer's warranty
- Ships with wheels installed and ready to use
- Minimal assembly required at delivery
- Latex free
- Durable, hard-baked powder coat finish

EMERGENCY Package Accessories Include:

- Cardiac board and brackets
- Powder coated steel top rail / push handle
- Side mounted drop shelf
- IV pole
- Oxygen tank holder (D/E size cylinders)

Internal Drawer Dimensions:

- 23" W x 16.75" D (58.4 cm W x 42.5 cm D)

Assembled Cart Dimensions

- 67.4" H x 33.25" W x 22" D (171.2 cm H x 84.5 W x 55.9 D) (height includes IV Pole)
- Cart weight: 138 lbs (62.6 kg)
- Approx shipping weight: 157 pounds (71.2 kg)