



AUTOMATIC Blood Pressure Monitor

UA-621

- Irregular HeartBeat (IHB) Indicator
- Blood Pressure Classification Indicator
- SlimFit Comfortable Cuff
- Fully Automatic Measurement
- Large 30-Memory + Average Reading
- Standard AC adapter Inlet

* Asia/EC model



A&D
A&D Medical
Technology You Can Trust

<http://www.aandd.jp>



WRIST Blood Pressure Monitor

UB-510

- Irregular HeartBeat (IHB) Indicator
- Blood Pressure Classification Indicator
- Comfortable Wrist Cuff
- Fast Measurement along with Inflation
- Large 30-Memory + Average Reading

* Asia/EC model



A&D
A&D Medical
Technology You Can Trust

<http://www.aandd.jp>

UA-621 Upper Arm Automatic Blood Pressure Monitor

UB-510 Wrist Blood Pressure Monitor

*Asia/EC model



Irregular HeartBeat Indicator

The indicator illuminates if an IHB is detected during the measurement. The indicator alerts users of an irregular heartbeat, which may be associated with cardiovascular disease. (Note: This is for warning purposes only. A&D recommends contacting a physician if the user observes this indicator frequently.)

Blood Pressure Classification Indicator

The World Health Organization (WHO) has designated 6 blood pressure classifications. This monitor automatically indicates your classification (Optimal to Severe hypertension) after each measurement is completed.

| Systolic Pressure | Diastolic Pressure | |
|-------------------|--------------------|---------------------------------|
| 180~ | 110~ | Grade 3 hypertension (severe) |
| 160~179 | 100~109 | Grade 2 hypertension (moderate) |
| 140~159 | 90~99 | Grade 1 hypertension (mild) |
| 130~139 | 85~89 | High-normal |
| 120~129 | 80~84 | Normal |
| Less than 120 | Less than 80 | Optimal |

Reference: Journal of Hypertension 1999, Vol 17 No.2

Multilingual Operation manual

| Asia model | | EC model | | |
|------------|------------|----------|------------|---------|
| English | Chinese | English | Malayalam | Arabic |
| Korean | Thai | Spanish | French | Chinese |
| Malay | Vietnamese | Italian | Portuguese | |
| Hindi | Indonesian | Turkey | German | |

SlimFit Comfortable Cuff (UA-621)

During a measurement, our new innovative edge free internal bladder "EFB" distributes pressure evenly around the user's arm. This provides a higher level of comfort while the arm is under pressure.

Fully Automatic Measurement

Pressing the START button is the only step necessary to perform a complete and accurate measurement.

Large 30-Memory + Average Reading

Both the UA-621 and the UB-510 have a 30-reading memory capacity for monitoring your measurement trends. An average of all the stored data is shown at the start of each measurement, offering the user their comprehensive measurement history at a glance.

Standard AC adapter Inlet (UA-621)

The monitor is equipped with an AC inlet for use with an external power source. Since there is no need to worry about battery life, the user can economize while enjoying the use of the monitor everyday.

(The AC adapter is optional and should be bought separately.)

A&D
A&D Medical
Technology You Can Trust
<http://www.aandd.jp>

Exclusively Distributed by:

WellnessPRO
 tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor. Lourdes St., Ortigas Ctr., Pasig City
 Tel: (02) 636-3590
 Fax: (02) 636-7442
 Mobile: (63) 920-947-7185, (63) 922-823-2282
 Email: wellnessproinc@gmail.com
www.wellnessproinc.com