




# FIBRE OPTIC (F.O.) LARYNGOSCOPES



**KaWe MEGALIGHT®  
MACINTOSH F.O.**




## Key Features:

- Matt finish stainless steel
- With integrated fibre optics
- No dust collecting cracks (easy to clean)
- Very large fibre optics
- Meaningful brightness
- Autoclavable
- Fits with KaWe and other DIN EN ISO 7376 handles
- Maintenance-free



## Other Blade Types:



*\*conventional / non-fibre optic models also available, call for more details"*

*\*see back for available sizes*

# OPTIONS

## Bulb Type

XENON  
LED

\*XENON (Standard for KaWe Blades)

## Power Source Options

Plain Battery  
Rechargeable Battery

## Accessories


MedCharge® 4000

Case  
Bulb

## Light Source

At Blade Head (Conventional)  
At Handle Top (Fibre Optics - FO)

## Upclose view of lamp / bulb for each F.O. blade type



## Blade Type

NAME	Available Size	IMAGE
*KaWe MEGALIGHT® Macintosh F.O.	0 1 2 3 4	
*KaWe MEGALIGHT® Miller F.O.	00 0 1 2 3 4	
Macintosh F.O.	0 1 2 3 4 5	
Miller F.O.	00 0 1 2 3 4	
*TEPRO® Macintosh F.O.	2 3 4	
Foregger F.O.	0 1 2 3 4	
*KaWe FLAPLIGHT® F.O. MEGALIGHT	2 3 4	
Polio Macintosh F.O.	3 4	
KaWe FLAPLIGHT® F.O.	2 3 4	

\*Conventional / Non FO Blades come in same shape / size range as the FO series except for (\*).

## Handle Size



\*Conventional Handles come in same size / type range as FO Handles except for A and B

## Laryngoscope Set

### NEONATAL



Miller F.O.  
Blades included:  
00, 0, 1

### PEDIATRICS



Miller / Macintosh F.O.  
Blades included:  
0, 1, 2

### ADULT



KaWe MEGALIGHT®  
Macintosh F.O.  
Blades included:  
2, 3, 4

### ADULT ECONOMY



Economy  
Macintosh F.O.  
Blades included:  
2, 3, 4

### PEDIATRICS + ADULT



Macintosh F.O.  
Blades included:  
2, 3, 4  
Miller F.O.  
Blades included:  
0, 1, 2, 3

See Wellness PRO or Kawe website for details on handle power options, conventional type range, actual sizes of blades, available accessories, etc.

\*Please indicate if you need FO or Conventional Handle when ordering

Available through:



Exclusive Philippine Distributor:

**WellnessPRO**  
tomorrow's medical devices

U-808 Richmond Plaza, San Miguel Avenue cor. Lourdes Street  
Ortigas, Pasig 1605 Philippines  
Tel.: 02 636-3580 • Fax: 02 6367542 • Cel.: 639209477185  
Email: info@wellnessproinc.com • Url: www.wellnessproinc.com