



Mini Doppler ES-100VX



**POCKET DOPPLER
WITH GREAT
SOUNDS**

Mini Doppler ES-100VX



FOR VASCULAR AND FETAL

- **HAND-HELD DESIGN**
- **GREAT SENSITIVITY**
- **WIDE PROBE SELECTION**
For user flexibility. (2, 4, 5, 8 and 10 MHz)
- **CONVENIENT PROBE ACTIVATION BUTTON**
Turns ES-100VX ON and OFF
- **AUTOMATIC POWER "OFF"**

PROBE SELECTION

VASCULAR

STANDARD



FETAL



CLINICAL APPLICATIONS

- Peripheral vascular procedures
- Fetal heart rate
- Penile & Digit systolic pressures
- Venous compressions
- Blood Pressure segmental studies
- Flow detection in recovery room

SPECIFICATIONS

Power	DC9V, A 9 volt alkaline battery	Dimensions	75(W)x140(D)x25(H) mm
Battery life	Approx. 3 hours Automatic shut-off	Weight	Approx. 270 grams (including probe)
Multi frequency	2, 4, 5, 8, 10 MHz	Accessories (optional)	Carrying case, camera strap, gel, headset, battery
Probe button	Power ON/OFF	Headset	Speaker cut off
Speaker output	200 mW or more Outputs (3.5mm jack)		



Kawasaki, Japan
www.hadeco.co.jp/



Tel Nos.: +63 (2) 8636 3580 | +63 (2) 8637 0930
Mobiles: +63 917 851 3948 | +63 919 082 6498
Email: helpdesk@wellnessproinc.com

#56 San Rafael St., Brgy. Kapitolyo,
Pasig City 1603 Philippines

www.wellnessproinc.com