

# Otoscope/Ophthalmoscope

®



## KaWe PICCOLIGHT E56 / US Ophthalmoscope

01.23566.021

- Battery handle (AA)
- 2.5V halogen bulb
- With correction lens wheel in stages from +/- 20 dioptres
- Aspherical precision optics
- In cloth bag
- Batteries: 2 Mignon, 1.5V (type AA)
- 6 diaphragms

6 apertures



## KaWe PICCOLIGHT E56 / US Ophthalmoscope

01.13600.022

- F.0 LED illumination
- Illuminance approx. 30,000 Lux
- Lifetime bulb approx. 100,000 h
- Energy-saving
- Battery handle AA
- Disposable ear funnels 10 x 0 2.5 mm
- Disposable ear funnels 10 x 0 4.0 mm
- In cloth bag
- Battery: 2 Mignon (type AA)



Disposable Ear Funnels Included