



KaWe MASTERMED® C

PREMIUM BLOOD PRESSURE MEASURING DEVICES



The KaWe MASTERMED® C Premium blood pressure monitors provide clinicians the highest level of quality and value. This series is defined by its innovative and practical features, smooth surfaces for easy cleaning, as well as its particularly precise gauge. These blood pressure monitoring devices are available as wall-mounted, table-mounted and freestanding models. These units have a double-walled housing that ensures additional stability and protection and are therefore designed to withstand twice the pressure. Each model comes with a generously-sized cuff holder.

Features

- ▶ complies with ISO 81060-1 standards
- ▶ clear gauge with thin, red dial for exact measurement readings
- ▶ numbered instrument dial reduces the risk of parallax
- ▶ 150 mm dial diameter
- ▶ scaling: 0 - 300 mmHg
- ▶ beryllium-copper diaphragm maintains the calibration up to 600 mmHg
- ▶ swivel-mounted
- ▶ equipped with PU/nylon-Velcro cuff for adults
- ▶ plug-in adapter for quicker cuff exchange
- ▶ stays calibrated for 5 years, guaranteed
- ▶ 2-year guarantee on device when used properly



Beryllium-copper diaphragm



Cuff holder



Metal valve



Velcro cuff



WALL-MOUNTED MODEL

- swivelling mounting bracket
- for use with the wall-mounted KaWe MedCenter® 5000 work station
- screws and wall anchors included
- spiral hose with a range of 3 m

TABLE TOP

- integrated handle and bases that prevent scratches
- spiral hose with a range of 1.2 m



STAND MODEL

- stable 5-foot roller stand made of high-quality stainless steel and sturdy plastic juncture for the feet
- adjustable height: 85 to 130 cm
- quiet Ø 50 mm double castors: 2 locking castors, 2 safety castors and one antistatic castor
- spiral hose with a range of 3 m

