

PROFI KARDIOLOGIE | PROFI CARDIOLOGY

Das Profi-Kardiologie Stethoskop ist das ideale Stethoskop für den Kardiologen. Die Bruststücke der KaWe Profi-Kardiologie Stethoskope sind in massiver Edelstahlausführung lieferbar. Anhand der hochwertigen Verarbeitung und Technik wurden Qualitäts-Stethoskope mit hervorragender Akustik geschaffen. Dieses hochwertige Stethoskop verfügt über ein Monoschlauchsystem mit integriertem Ohrbügel aus Edelstahl, weiche Softoliven sorgen für einen optimalen Tragekomfort. Jedem Stethoskop liegen zwei Paar Ersatz-Ohrliven, eine Ersatz-Membrane und ein Namensschild bei.

The Profi Cardiology stethoscope is the ideal stethoscope for cardiologists. The chest piece of the KaWe Profi Cardiology stethoscope is available in solid stainless steel. These superior stethoscopes and their outstanding acoustics are crafted with high-quality workmanship and technology. This high-quality stethoscope has a single-tube system with an integrated head piece made of stainless steel. Soft ear tips ensure optimal wearer comfort. Each stethoscope comes with two spare pairs of ear tips, one spare diaphragm and a name plate.

Profi-Kardiologie Stethoskop Profi Cardiology stethoscope

- Doppelkopf-Stethoskop
 - Bruststück in massiver Edelstahlausführung
 - Membran Ø 45 mm
 - Schlauchlänge ca. 42 cm
 - umschaltbar
 - Monoschlauchsystem mit integriertem Ohrbügel aus Edelstahl
 - 2 Paar Ersatz-Ohrliven, 1 Ersatz-Membrane, 1 Namensschild
- dual-head stethoscope
 - chest piece in solid stainless steel design
 - diaphragm Ø 45 mm
 - tube length approx. 42 cm
 - reversible
 - single-tube system with integrated head piece of stainless steel
 - 2 pairs of spare ear tips, 1 spare diaphragm, 1 name plate

schwarz black	REF 06.22800.022 (43440)
burgunder burgundy	REF 06.22800.172 (43441)



138

ERSATZTEILE | SPARE PARTS

für Profi- und Top-Kardiologie | for Profi and Top Cardiology

Ohrliven weich, schraubbar
Ear tips soft, screw-on

schwarz, VE = 10 Stk.
black, PU = 10 items

REF 06.91220.021

Ohrliven weich, schraubbar
Ear tips soft, screw-on

burgunder, VE = 10 Stk.
burgundy, PU = 10 items

REF 06.91220.171

Ohrbügel mit Schlauch
Head piece with tube

Schlauchlänge ca. 42 cm
tube length approx. 42 cm

schwarz black	REF 06.94280.021 (42620)
burgunder burgundy	REF 06.94280.171 (42621)

**KARDIOLOGIE
CARDIOLOGY**

~35 mm

~45 mm