


WIN THE WAR ON CROSS INFECTION!



KaWe SINGLE PATIENT USE (SPU) CUFF

*Space for Hospital Logo - Minimum Quantity Required

A. INFANT SIZE KAWE-DC-DT-INF



○ arm size: 13.8 - 21.5cm

B. CHILD SIZE KAWE-DC-DT-CHI



○ arm size: 20.5 - 28.5 cm

C. ADULT SIZE KAWE-DC-DT-ADU



○ arm size: 27.5 - 36.5 cm

D. LARGE ADULT SIZE KAWE-DC-DT-LAD



○ arm size: 35.5 - 46.0 cm



Key Features:

- Suitable for Long-Term Use
- Patient Identification
- Prevents Cross-Contamination
- Comfortable Fleece Material
- White Color for Easy Inspection
- Size Guide

Key Benefits:

- Fits into Daily Nursing Care
- Saves Time
- Saves Money
- Good for All Manometers
- Comfortable and Safe for Patient

Exclusively Distributed by:


WellnessPRO
tomorrow's medical devices

Readily Fits **KaWe** Blood Pressure Heavy Duty Device

- Okay up to 600mgHg with no damage
- Copper Bellirium Internal Workings



KaWe MC100
Table Top




KaWe MC102
Wall Mounted



KaWe MC101
Wheeled Stand

It's all about Easy-Connect



Disposable Cuff with
Connectors at end of
Tube A and B



Rubber Bulb with
Built-In Valve Having
Connector

Infection Control Merits of Single Use **KaWe** Blood Pressure Cuff

1. Reusable BP Cuff has been found to be a major source of contamination, especially in Intensive Care Units leading to Hospital Acquired Infection (HAI).
2. WHO Guidelines require regular cleaning of Blood Pressure Cuffs.
3. Reduction in Hospital Acquired Infection (HAI) cases directly impacts your hospital's staff workload as well as the overall cost of operations.

Technical Specifications of **KaWe** Cuffs

Single Tube with Connector Tube

Sequence:	1	2	3	4
Item #:	KaWE DCCT-ST-INF	KaWE DCCT-ST-CHI	KaWE DCCT-ST-ADU	KaWE DCCT-ST-LAD
Description:	Kawe Disposable Cuff, Single Tube, Infant, with connector tube	Kawe Disposable Cuff, Single Tube, Child, with connector tube	Kawe Disposable Cuff, Single Tube, Adult, with connector tube	Kawe Disposable Cuff, Single Tube, Large Adult, with connector tube
Dimensions:	13.8-21.5cm	20.5-28.5cm	27.5-36.5cm	35.5-46.0cm
Packing:	10	10	10	10
Photo:				

Double Tube with Connector Tube

Sequence:	1	2	3	4
Item #:	KaWE DCCT-DT-INF	KaWE DCCT-DT-CHI	KaWE DCCT-DT-ADU	KaWE DCCT-DT-LAD
Description:	Kawe Disposable Cuff, Double Tube, Infant, with connector tube	Kawe Disposable Cuff, Double Tube, Child, with connector tube	Kawe Disposable Cuff, Double Tube, Adult, with connector tube	Kawe Disposable Cuff, Double Tube, Large Adult, with connector tube
Dimensions:	13.8-21.5cm	20.5-28.5cm	27.5-36.5cm	35.5-46.0cm
Packing:	10	10	10	10
Photo:				

Single Tube without Connector Tube

Sequence:	1	2	3	4
Item #:	KaWE DC-ST-INF	KaWE DC-ST-CHI	KaWE DC-ST-ADU	KaWE DC-ST-LAD
Description:	Kawe Disposable Cuff, Single Tube, Infant, without connector tube	Kawe Disposable Cuff, Single Tube, Child, without connector tube	Kawe Disposable Cuff, Single Tube, Adult, without connector tube	Kawe Disposable Cuff, Single Tube, Large Adult, without connector tube
Dimensions:	13.8-21.5cm	20.5-28.5cm	27.5-36.5cm	35.5-46.0cm
Packing:	10	10	10	10
Photo:				

Double Tube without Connector Tube

Sequence:	1	2	3	4
Item #:	KaWE DC-DT-INF	KaWE DC-DT-CHI	KaWE DC-DT-ADU	KaWE DC-DT-LAD
Description:	Kawe Disposable Cuff, Double Tube, Infant, without connector tube	Kawe Disposable Cuff, Double Tube, Child, without connector tube	Kawe Disposable Cuff, Double Tube, Adult, without connector tube	Kawe Disposable Cuff, Double Tube, Large Adult, without connector tube
Dimensions:	13.8-21.5cm	20.5-28.5cm	27.5-36.5cm	35.5-46.0cm
Packing:	10	10	10	10
Photo:				

Product Specifications and Technical Details may change.

Exclusively Distributed by:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Avenue corner Lourdes St., Ortigas Center Pasig 1605, Philippines

Tel: (02) 636-3580 • Fax: (02) 636-7542
Mobile: (63) 920-947-7185, (63) 922-823-2282

Email: wellnessproinc@gmail.com
jcjuan@wellnessproinc.com
www.wellnessproinc.com