



## RX Series

Hand hygiene made easy

Exclusively Distributed by:

**WellnessPRO**  
tomorrow's medical devices

## RX Series - Hand hygiene made easy

Manual or touchless – the new RX dispenser line is a flexible and economical range for any medical setting. Whether you're looking for a dispenser for intensive care units or the hospital lobby, the RX series offers wall mounting, countertop/trolley (only RX 5 T), or hygiene stand solutions - enabling proper hand hygiene in every location. With a robust plastic housing, **compatibility with disposable hygiene pumps and packaging (conforming with the latest hygiene regulations and standards)**, and optional hygiene station accessories, this dispenser series offers everything you need for a complete hygiene solution.



Soap and disinfectant dispenser, 500 ml

### RX 5 M

for disinfectant, alco-gel, soap and flowable lotions, for 500 ml Euro bottle, accepts disposable hygiene pumps and packaging, plastic, incl. exchangeable pump 1.5 ml / stroke, stainless steel lever, drip tray (removable), for wall mounting, incl. HDPE bottle (with standard pump option), with optional lock.

Standard pump:	Dosage amount 1.5 ml / stroke
Disposable hygiene pump:	Dosage amount 1.5 ml / stroke
Bottle:	For 500 ml Standard Euro bottle, refillable or replaceable with any product adhering to Euro Bottle standard
Dimensions:	W 98 (176) x H 300 x D 144 (192) mm (with drip tray)

	with Standard Pump	Art. No.
RX 5 M (Non-Locking)	with bottle	4401841
RX 5 M (Locking)	Lockable, with bottle	4401842
	with Disposable Hygiene Pump	Art.-No.
RX 5 M DHP (Non-Locking)	without bottle	4401974
RX 5 M DHP (Locking)	Lockable, without bottle	4401975

„Safety-Click“ Drip tray

Exclusively Distributed by:

**WellnessPRO**  
tomorrow's medical devices



Touchless soap and disinfectant dispenser, 500 ml

**RX 5 T**


for disinfectant, alco-gel, soap and flowable lotions, for 500 ml Euro bottle, accepts disposable hygiene pumps and packaging, plastic, incl. exchangeable pump 1.0 or 1.5 ml / stroke, mobile device, optional with back plate for wall mounting, hygiene stand ordered separately incl. HDPE bottle and 4 C cell batteries.

- Pump: Dosage amount 1.0 or 1.5 ml / stroke, can be varied by selecting one, two or three strokes per activation
- Bottle: For 500 ml Standard Euro bottle, refillable or replaceable with any product adhering to Euro Bottle standard
- Dimensions: W 102 x H 269 x D 142 (158) mm (with back plate)

	<b>1.0 ml / stroke dosage</b>	<b>Art. No.</b>
RX 5 T	Tabletop device	4401736
RX 5 T	Tabletop device with back plate for wall mounting	4401764
	<b>1.5 ml / stroke dosage</b>	<b>Art. No.</b>
RX 5 T	Tabletop device	4401117
RX 5 T	Tabletop device with back plate for wall mounting	4401528
	<b>Optional Accessories</b>	<b>Art. No.</b>
Back plate	For wall mounting	4401226



Exclusively Distributed by:



Costly reprocessing eliminated!

### RX 5 Pumps

#### DHP

Disposable hygiene pump for manual and touchless RX dispensers

Single-use plastic pump reduces the need for cleaning/sterilization reprocessing.

Material chemically compatible with a diverse range of liquids.

Minimum order quantity one box (150pcs.).

Dosage: 1.5 ml/ stroke



	Disposable Hygiene Pump	Art.-No.
▶	DHP RX	4401860



#### IMP K

Standard plastic pump for manual and touchless RX dispensers


High-quality, long-life plastic pump.

For wide range of soaps, lotions and disinfectants, material chemically compatible with a diverse range of liquids, dishwasher safe.

Dosage: 0.7 / 1.0 / 1.2 / 1.5 ml/ stroke

Cleaning: Dishwasher safe (60°C)

	Standard Plastic Pump	Art.-No.
▶	IMP K for RX 5 M/T	4401434



### Stainless steel stands

#### RX Stands

Robust stainless steel stands

Offers convenient hand hygiene locations in nearly every situation

	for RX 5 M	Art.-No.	
▶	Hygiene Station E AFP	Station without dispenser and drip tray	1420177
	Document holder	DIN A4, retrofittable	1420527
	for RX 5 T	Art.-No.	
▶	RX 5 T Stand	Height 122 cm, stand ø 40 c, with lock, without dispenser	1421065

Exclusively Distributed by: