

# seca 700

## Mechanical column scales with eye-level beam



- Extra high capacity of 220 kg, also for heavyweight patients.
- Large, low platform.
- Simple to operate and easy to read due to eye-level display.
- Sturdy construction with cast-iron base.
- Transport castors for mobile use.

## seca 700:

### Precision and quality with character.

**A**bsolute precision and an indestructible mechanism are prerequisites for a scale which is used daily in professional health care. With these qualities, the reliable column scale with the characteristic beam at eye level has made a name for itself in medical practice. The optimized, updated model seca 700 now has an even higher capacity of 220 kg, is easier to operate and has an improved design. Its profile is in keeping with the modern requirements in general medicine, in internal medicine and in nutrition consultation.



#### High capacity and large platform.

With its capacity of 220 kg, the seca 700 has been specially designed for weighing overweight and adipose persons. To ensure durability, special attention was paid to the high-quality construction – right down to the smallest detail. The generously sized, low level platform make the seca 700 easy to mount. Weighing has thus become much more comfortable, not only for the patients but also for doctors and carers.

#### Easy to operate and read due to eye-level display.

The weight scale, which has been affixed at eye level so that the result can be read off comfortably, is printed on both sides so that it can be read both by the person standing on the scale and by the person standing in front of the scale. The large numerals make it easier to quickly determine the weight. A special coating reliably prevents the scale from fading even after years of use.

#### Sturdy construction of cast iron and integrated transport castors.

The base of the seca 700 is made of cast iron and can cope with even the most intensive usage. The special coating on the lever apparatus protects the inner parts against corrosion and therefore guarantees precise results for a long time. In spite of its compact construction, the scale can be moved easily and quickly on the integrated transport castors to wherever it is needed. Fitted with the optional telescopic measuring rod seca 220, the seca 700 makes time-saving weighing and measuring possible in one step.



The optional measuring rod seca 220 allows convenient weighing and measuring in one step.



The broad and flat construction provides easy step-on access and a firm foothold.

## seca 700

### Technical Data

- Capacity: 220 kg, 500 lbs or 160 kg / 400 lbs
- Graduation: 50 g, 1/8 lbs or 100 g / 1/4 lbs
- Dimensions (WxHxD): 520 x 1,556 x 520 mm / 20.5 x 61.3 x 20.5"
- Dimensions platform (WxHxD): 335 x 80 x 345 mm / 13.2 x 3.1 x 13.6"
- Weight: 16.2 kg / 35.7 lbs
- Functions: reset-to-zero function,
- Optional: measuring rod seca 220

Exclusive Philippine Distributor:

**WellnessPRO**  
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,  
Ortigas Ctr. Pasig City 1605 Philippines  
Mobile no.: +63 922 896 0964  
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473  
Fax: +63 2 636 7542  
Email: wellnessproinc@gmail.com, icjuan@wellnessproinc.com  
[www.wellnessproinc.com](http://www.wellnessproinc.com)