

seca 760

Mechanical personal scales in modern colours



- Robust construction.
- Large, easy-to-read dial.
- Low profile and wide platform.
- Precise and durable.

seca 760:

A classic in several colours.

This seca classic has greatly influenced the image of scales. With an absolutely precise dial, synthetic covering used in medical applications, and steel plate housing from the good old days. The modern age has also left its mark on these scales, which are available in many different colours and models. For one thing a classic should never be boring.

In several colours to suit every taste.

These scales are available in burgundy red, midnight blue and ecru. And because there's no disputing taste, this classic model also comes in two additional versions – with white or black housing. All of them have

in common the chrome-plated safety ring, a platform of skin-friendly grooved natural rubber, the classic scale dial and robust, precise weighing mechanics.



Traditional weighing mechanics in classic colours.



Exclusive Philippine Distributor:

WellnessPRO
tomorrow's medical devices

Unit 808 Richmond Plaza, San Miguel Ave cor Lourdes St.,
Ortigas Ctr. Pasig City 1605 Philippines
Mobile no.: +63 922 896 0964
Tel. no.: +63 2 636 3580, +63 2 637 0930, +63 2 637 0473
Fax: +63 2 636 7542
Email: wellnessproinc@gmail.com, jjcjuan@wellnessproinc.com
www.wellnessproinc.com

seca 760

Technical Data

- Capacity: 150 kg, 320 lbs
- Graduation: 1.000 g, 1 lbs
- Dimensions (WxHxD):
303 x 118 x 470 mm
- Weight: 3.5 kg / 7.7 lbs
- Functions: reset-to-zero function